


"Collect & Sterilization of N95 Respirator Masks":

<u>Purpose:</u> To establish a process to collect and sterilize N95 Respirator Masks across the system to ensure we have

adequate supply.

Scope: Applicable in all care settings.

Approved: 4/9/2020 Revised: 4/9/2020



SPD will process N95

mask and specify the

frequency as which the

units come to SPD

The designated person, on each unit, will go to

SPD to pick up the processed masks for distribution to the staff

## NOTE:

- 1. DO NOT USE: Makeup, Lipstick, lip balm, Chap Stick, Vaseline, facial cream or lotion.
- 2. If the masks is wet or soiled then it cannot be reprocessed.
- 3. Per FDA regulations, the masks will only be able to be processed 2 times.
- 4. Reprocessing- The systematic approach used to ensure the pathogens are killed on the mask including COVID-19 following guidelines set forth by the FDA.
- 5. You will get your own mask returned to you.